# **ECG** Interpretation using Python

1. Importing important libraries:

```
import os
import pandas as pd
import pandas as pd
import pandas as pd
import numpy as np
import matplotlib.pyplot as plt
import math
import os
import re
import scipy
from scipy.signal import find_peaks, find_peaks_cwt
```

## **Reading Data**


2. Reading the data from the CSV file: This will read the file and return a dataframe. And it also makes sure that file is in CSV or XLSX format.

```
In [ ]: def read file(file path):
            Read a file either as CSV or Excel based on the file extension.
            Parameters:
            - file_path (str): Path to the file.
            - pd.DataFrame: DataFrame containing the data from the file.
            if file_path.endswith('.csv'):
                df = pd.read_csv(file_path)
            elif file_path.endswith('.xlsx') or file_path.endswith('.xls'):
                df = pd.read_excel(file_path)
            else:
                return None
            return df
        file path = r'C:\Users\Divyansh\Desktop\U4RAD\READ MY ECG\ALPHA 2023.xlsx'
        DF = read_file(file_path)
        DF.columns = DF.columns.str.strip()
        column_name = 'II'
```

3. Selecting the data from Lead II and plotting the graph.

```
In [ ]: ekg = DF[column_name].copy()
    ekg.plot()
    plt.axhline(y=0, color='r', linestyle='--', label="y=0")
```


```
plt.title("EKG")
plt.show()
```



# Filtering the Data


4. Filtering the data and smoothening it using the Savgol Filter from scipy.

```
In [ ]: smoothed_heartbeats = scipy.signal.savgol_filter(ekg, window_length=20, polyorde
    smoothed_heartbeats = pd.Series(smoothed_heartbeats)
# print(smoothed_heartbeats)
smoothed_heartbeats.plot()
plt.axhline(y=0, color='r', linestyle='--', label="y=0")
plt.title("Smoothed_EKG")
plt.show()
```



5. Flattening the graph as to correct the wandering of the graph using a medium filter medfilt() from scipy.

```
In []: frequency = 500 # Hz,
    kernel_size = frequency + 1
    wandering_baseline = scipy.signal.medfilt(ekg, kernel_size=kernel_size)
    flattened_ekg = ekg - wandering_baseline
    flattened_ekg.plot()
    plt.axhline(y=0, color='r', linestyle='--', label="y=0")
    plt.title("Flattened EKG")
    plt.show()
```



#### **RR Interval and Heart Beat**

6. Calculating RR Intervals and Heart Beat per minute from the data.

Function to calculate heart beat per minute by taking RR Interval as an arguement.

```
if(interval==0):
                hbpm = "0"
                return hbpm;
            number_of_blocks = interval/200
            hbpm = 300/number_of_blocks
            return math.floor(hbpm)
In [ ]: df = pd.DataFrame({'ecg': flattened_ekg})
        column_name = 'ecg'
        # Find the R-points.
        # Find peaks in the signal
        peaks, _ = find_peaks(df[column_name])
        # Find troughs in the signal (minima)
        troughs, _ = find_peaks(-df[column_name])
        try:
              positive_peaks = [peak for peak in peaks if df[column_name][peak] > 0]
              # Filter peaks with values greater than 0.2
              selected_peaks = [peak for peak in peaks if df[column_name][peak] > 0.35]
              # Calculate the average interval between maxima
              if len(selected_peaks) > 1:
                 average_interval = sum(selected_peaks[i+1] - selected_peaks[i] for i in
```

else:

def heart\_beat(interval):

In [ ]:

```
average_interval = 0

# Calculate the average value of the maxima
average_maxima_value = sum(df[column_name][peak] for peak in selected_peak

R_II = "{:.2f}".format(average_maxima_value)
RR_Interval=math.floor(average_interval)*2
hbpm = heart_beat(RR_Interval)


except Exception as e:
    print(e)

print("RR Interval:",RR_Interval)
print("R(II):",R_II)
print("HR:",hbpm)
```

RR Interval: 770 R(II): 1.07 HR: 77

7. Plotting the R-peaks in the graphs.

```
In []: plt.plot(df[column_name])
    plt.plot(selected_peaks, df[column_name][selected_peaks], "x", label="R Points",
    plt.axhline(y=0, color='r', linestyle='--', label="y=0")
# Adding labels and title
    plt.xlabel('X-axis label')
    plt.ylabel('Y-axis label')
    plt.title('Line Chart for {}'.format(column_name))
# Display the chart
    plt.legend()
    plt.show()
```



## **QRS Complex**

8. Recognizing the QRS complex and finding out QRS duration.

```
In [ ]: #QRS
        try:
            # Initialize lists to store Q and S points
            q_points = []
            s points = []
            # Identify Q and S points for each R point
            for r peak in selected peaks:
                 # Find troughs before and after the R peak
                troughs_before_r = [trough for trough in troughs if trough < r_peak]</pre>
                troughs after r = [trough for trough in troughs if trough > r peak]
                 # Find the Q point (minima with negative values just before R)
                q_point = max((trough for trough in troughs_before_r if df[column_name][
                # Find the S point (minima with negative values just after R)
                 s_point = min((trough for trough in troughs_after_r if df[column_name][t
                 # Append Q and S points to the respective lists
                 q_points.append(q_point)
                 s_points.append(s_point)
            newQ_points = []
            for point in q points:
                 # Find the peak just to the left of the point
                 newQ_point = max((peak for peak in peaks if peak < point), default=None)</pre>
```

```
# Append the found peak to the list
        newQ_points.append(newQ_point)
    newS_points = []
    for point in s_points:
        # Find the peak just to the left of the point
        newS_point = min((peak for peak in peaks if peak > point), default=None
        # Append the found peak to the list
        newS_points.append(newS_point)
    QRS durations = []
    for q_point, s_point in zip(newQ_points, newS_points):
        if q_point is not None and s_point is not None:
            duration = s_point - q_point
            QRS_durations.append(duration)
    # Calculate average duration
    average_duration_qrs = sum(QRS_durations) / len(QRS_durations) if len(QRS_du
    QRS_avg = math.floor(average_duration_qrs*2)
except Exception as e:
   print(e)
   QRS_avg=0
    # print("QRS BLOCK",patient_name,e)
print("QRS:",QRS_avg)
```


QRS: 81

9. Plotting the QRS Complex.

```
In []: plt.plot(df[column_name])
    plt.plot(selected_peaks, df[column_name][selected_peaks], "x", label="R Points",
    plt.plot(newQ_points, df[column_name][newQ_points], "o", label="Q Point", color=
    plt.plot(newS_points, df[column_name][newS_points], "o", label="S Point", color=
    plt.axhline(y=0, color='r', linestyle='--', label="y=0")

# Adding LabeLs and title
    plt.xlabel('X-axis label')
    plt.ylabel('Y-axis label')
    plt.title('Line Chart for {}'.format(column_name))

# Display the chart
    plt.legend()
    plt.show()
```



# **QT Interval**

10. Finding T-peaks and calculating values of QT, QTc and QT/QTc Ratio.

```
In []: def find_maximum_between_ranges(r_point, p_point):
    max_values = []

for i in range(len(r_point) - 1):
    start_index = p_point.index(r_point[i])
    end_index = p_point.index(r_point[i + 1])
    values_within_range = p_point[start_index + 1:end_index]

    if values_within_range:
        max_value = max(values_within_range)
        max_values.append(max_value)
    else:
        max_values.append(None)

return max_values
```

```
In []: r_points = []
    for r_point in selected_peaks:
        r_points.append(df[column_name][r_point])
    r_pair = [(key,value) for i, (key,value) in enumerate(zip(selected_peaks, r_poin r_dict = dict(r_pair)
    #print(r_dict)

p_points = []
    for p_point in positive_peaks:
```

```
p_points.append(df[column_name][p_point])
# print("P-points",p_points)
p_pair = [(key,value) for i, (key,value) in enumerate(zip(positive_peaks, p_point p_dict = dict(p_pair))
# print("P DICT:",p_dict)
```

```
In [ ]: t_point = find_maximum_between_ranges(r_points,p_points)
        # Create a dictionary where keys are taken from p_dict and values from t_point
        t_dict = {key: value for key, value in p_dict.items() if p_dict[key] in t_point}
        #print("T Dictionary:", t_dict)
        t_keys_list = [key for key in t_dict]
        #print("T Keys List:", t_keys_list)
        wave_end=[]
        for key in t_keys_list:
            for i in range(key,len(df.values)):
                if (df[column_name][i]<0):</pre>
                     wave_end.append(i)
                     break;
            # qt_durations = [t - q for q, t in zip(q_points, wave_end)]
        QT_durations = []
        for q_point, t_point in zip(newQ_points, wave_end):
            if q_point is not None and t_point is not None:
                duration = t_point - q_point
                 QT_durations.append(duration)
        avg_duration = (sum(QT_durations))/len(QT_durations))*2
        QT_avg = math.floor(avg_duration)
        QT_max = max(QT_durations)*2
        QT_min = min(QT_durations)*2
        if QT_avg>500:
            QT_avg = QT_min
        # print(QT_durations)
        # QT_interval = max(qt_durations)*2
        QTC = math.ceil(QT_avg/ math.sqrt(60/hbpm))
        QT_QC_Ratio = "{:.2f}".format(QT_avg/QTC)
        print("QT:",QT_avg)
        print("QTc:",QTC)
        print("QT/QTc:",QT_QC_Ratio)
```

QT: 347 QTc: 394 QT/QTc: 0.88


11. Plotting the wave end of the heartbeat.

```
In []: plt.plot(df[column_name])
    plt.plot(selected_peaks, df[column_name][selected_peaks], "x", label="R Points",
    plt.plot(newQ_points, df[column_name][newQ_points], "o", label="Q Point", color=
    plt.plot(newS_points, df[column_name][newS_points], "o", label="S Point", color=
    plt.plot(wave_end, df[column_name][wave_end], "o", label="End", color='black')
    plt.axhline(y=0, color='r', linestyle='--', label="y=0")

# Adding LabeLs and title
    plt.xlabel('X-axis label')
    plt.ylabel('Y-axis label')
    plt.title('Line Chart for {}'.format(column_name))

# Display the chart
```

```
plt.legend()
plt.show()
```



#### **PR Interval**

12. Finding the P-points and Calculating PR Interval.

```
In [ ]:
        def find_max_peak_in_range(df, column_name, start_index, end_index):
                # Ensure that the range is valid
                if start_index < 0 or end_index >= len(df):
                    print("Invalid range.")
                    return None
                # Extract the specified range from the DataFrame column
                range_data = df[column_name][start_index:end_index+1]
                # Find peaks within the range
                peaks = [(peak, start_index + i) for i, peak in enumerate(range_data)]
                if not peaks:
                    return None # No peaks found, return None
                else:
                    # Find the maximum peak within the range
                    max_peak = max(peaks, key=lambda x: x[0])
                return max peak
In [ ]:
         #Calculation of NEW PR
            # print("Wave-End", wave end)
            # print("Q-points", newQ_points, len(newQ_points))
        newQ_points.pop(0)
```

```
# print("newQ",newQ points)
P_points = []
for i in range(len(wave_end)):
    point = find_max_peak_in_range(df, column_name, wave_end[i], newQ_points[i])
    P points.append(point)
# print("P-points",P_points)
# print("P-points",P points)
P_pair = dict(filter(lambda x: x is not None, P_points))
# print("P-Pair:",P_pair)
values_list = list(P_pair.values())
    # print("P-peaks", values list)
p_start=[]
for p_peak in values_list:
   troughs_before_p = [trough for trough in troughs if trough < p_peak]</pre>
        # Find the Q point (minima with negative values just before R)
   start_point = max((trough for trough in troughs_before_p if df[column_name][
   p_start.append(start_point)
   # print(p_start)
    # P_start and new Q list should be equal in length
P_dict = {i + 2: value for i, value in enumerate(p_start)}
# print("P-DICT:",P_dict)
if len(newQ_points)>len(p_start):
   newQ_points.pop(0)
PR_durations = []
# result_dict = {key:Q_dict.get(key, None) - P_dict.get(key, None) if P_dict.ge
# print("P-Q:",result_dict)
# result_values = [value for value in result_dict.values() if value is not None]
# print("PR:",result values)
for p_point,q_point in zip(p_start, newQ_points):
    if q_point is not None and p_point is not None:
        duration = q_point - p_point
        PR_durations.append(duration)
PR avg= math.floor(sum(PR durations))/len(PR durations))*2
print("PR:",PR_avg)
```


PR: 132

#### 13. Plotting all the points.

```
In []: plt.plot(df[column_name])
    plt.plot(selected_peaks, df[column_name][selected_peaks], "x", label="R Points",
    plt.plot(newQ_points, df[column_name][newQ_points], "o", label="Q Point", color=
    plt.plot(newS_points, df[column_name][newS_points], "o", label="S Point", color=
    plt.plot(wave_end, df[column_name][wave_end], "o", label="End", color='black')
    plt.plot(p_start, df[column_name][p_start], "o", label="P-Start", color='green')
    plt.axhline(y=0, color='r', linestyle='--', label="y=0")

# Adding Labels and title
    plt.xlabel('X-axis label')
    plt.ylabel('Y-axis label')
    plt.title('Line Chart for {}'.format(column_name))

# Display the chart
    plt.legend()
    plt.show()
```



# **Making Data Set**

14. Making a Dictionary of all the data that we calculated.

```
In [ ]:
        data = {
                 "HR": hbpm,
                 "R(II)": R_II,
                 "RR": RR_Interval,
                 "PR": PR_avg,
                 "QRS": QRS_avg,
                 "QT": QT_avg,
                 "QTC": QTC,
                 "QT/QC": QT_QC_Ratio
        new_df = pd.DataFrame(columns=['HR','R(II)','RR','PR','QRS','QT','QTC','QT/QTC']
        row = list(data.values())
        new_df.loc[len(new_df)] = row
        print(new_df)
          HR R(II)
                          PR
                              QRS
                                     QΤ
                                         QTC QT/QTC
       0 77 1.07 770 132
                               81
                                   347
                                         394
                                               0.88
```